

# Welcome to our Quick Guide to Wearing the MOLLii Suit



*Always read the Instructions for Use*

## Putting a MOLLii Suit On



This presentation summarises how to put on a MOLLii Suit.


You must also read the *Instructions for Use* supplied with your MOLLii Suit.

## What you should have...



The Jacket, Trousers, Belt and Magnetic Clip Buckle and two Laundry Bags.

PLUS



The Programmed Control Unit with backplate and 4 x AAA batteries.

## Also...



A cross-head screwdriver, such as a No.1 Phillips Head or Pozidrive to screw the backplate onto the Control Unit.

## Start with putting the Trousers on...



- Put the stiffer or more sensitive leg in first.
- Ensure that the lower zip ends up on the outside of the leg and forms a line between the ankle and the outside of the knee.
- Make sure any socks are outside the trousers.

**Tip** : Carefully scrunch up a trouser leg to help it on especially for children.

Always read the *Instructions for Use*

## The Trousers continued.....

- Check that the zips on the trousers form a straight line on the outside of the leg.
- Check that the curved seam is above the kneecap.



## What it should look like...



- Pull the trousers right up before you do up the side zips and the fly.
- Ensure nappies, diapers or larger underwear do not obstruct the electrodes contacting the skin inside.

Socks outside trousers.



## Now put the Jacket on...



- Put the jacket on by inserting one arm at a time.
- Put the stiffer or more sensitive arm in first.
- Ensure that the side zips end up on the outside of the arm.

Always read the *Instructions for Use*

## Checking the Jacket...



- Check that the zips on the arms are positioned as in the picture.
- The elbow should be between the arm zips.
- The shirt should be outside the trousers.
- Do up the zip on the stomach first, then the zips for the arms.

## The Control Unit 1



- Insert the batteries.
- Thread the belt through the loops on the back of the control unit.
- Replace the backplate and secure the control unit and belt around the waist.

*Always read the Instructions for Use*

## The Control Unit 2




- Put the magnetic couplings on the contact pads on the garment.
- Ensure that they are positioned correctly.

## The Control Unit 3



Ensure the magnetic couplings are lined up with the contact pads on the garment and don't slip in use.



## Starting the Stimulation Session



Always read the *Instructions for Use*


- Ensure that the garment is dry.
- Start the control unit using the ON/OFF button:
- Press PLAY/PAUSE
- Once :



# Starting the Stimulation Session



The display shows:



The figures on the display show the remaining time in minutes.


## Starting the Stimulation Session



IF NECESSARY:  
PAUSE by pressing  
PLAY/PAUSE once :



The display then shows:



Press PLAY/PAUSE again  
to continue.

## At The End of the Session



- **FINISH:** Turn off the control unit using the ON/OFF button : 
- Remove the control unit before taking off the jacket and trousers.

*Always read the Instructions for Use*



For Further Enquiries:

Mollii Australia Head Office

416 High Street Maitland NSW 2320

PO Box 734 Maitland NSW 2320

P: 1300 315530

M: +61 491 616 452

E: [info@molliiaustralia.com.au](mailto:info@molliiaustralia.com.au)

W: [www.molliiaustralia.com.au](http://www.molliiaustralia.com.au)